

S/C.NEST™ | Investigation of seeding density and mixing speed for optimal CHO cell culture in 24-well scale mixing culture

Jess Lin, BS, James Chao, MS, Shih-Pei (Betty) Lin, PhD
CYTENA Bioprocess Solutions, Taipei City, Taiwan



Abstract

C.NEST and S.NEST are 24-well scale microbioreactors that offer the patented reciprocal mixing conditions for cell cultures. To obtain optimal mammalian cell growth in C/S.NEST, new users need to determine appropriate initial seeding density and mixing speed. In this study, we evaluated six initial seeding densities and three mixing speeds for CHO cell 24-well culture in C/S.NEST. Our findings offer guidance for researchers looking to incorporate the NEST system into their cell culture workflows.

Introduction

Chinese Hamster Ovary (CHO) cells are the workhorse of the biopharmaceutical industry, playing a pivotal role in the production of therapeutic proteins and monoclonal antibodies [1]. Over the past few decades, as the demands for biopharmaceuticals have surged, there has been a concomitant emphasis on refining cell culture methodologies to achieve maximal product yields [2]. One of the foremost parameters under investigation in this quest for optimization is seeding density.

Seeding density sets the trajectory for nutrient consumption, metabolic waste accumulation, and overall cell health in the bioreactor environment [3]. An improper seeding density can lead to reduced growth rates, sub-optimal protein production, and even apoptosis. Furthermore, studies have demonstrated that cell-seeding density is linked to the production yield of monoclonal antibodies in CHO cells, underscoring the necessity to control this parameter [4].

In this context, we performed this study to explore how initial seeding densities and mixing culture speeds may influence CHO-S cell growth in a NEST microbioreactor, a system that replicates bioreactor conditions at a compact 24-well scale.

Materials and methods

Cell culture

A monoclonal antibody (mAb)-expressing CHO-S cell line was used in this study. CHO-S cells were cultured in CD Hybridoma Medium (Thermo Fisher Scientific, Gibco, #11279023) supplemented with 100 units/mL penicillin, 100 µg/mL streptomycin (Thermo Fisher Scientific, Gibco, #5140122), 8 mM L-Glutamine (Corning, 25-005-CI), 0.2% Anti-Clumping Agent (Thermo Fisher Scientific, Gibco, #0010057DG) and 1X Cholesterol Lipid Concentrate (Thermo Fisher Scientific, Gibco, #12531018). The cells were cultured in either a humidified standard incubator or S.NEST with 5% CO₂ at 37°C.

Mixing culture with different cell seeding density

CHO-S cells were prepared in 2x10⁴, 5x10⁴, 1x10⁵, 2x10⁵, 5x10⁵ and 1x10⁶ cells/mL concentrations for seeding into four 24-well plates (Greiner, #662160). 1.4 mL of cell suspensions at each cell density was dispensed into 4 wells in each 24-well plate. 18 mL of Dulbecco's phosphate-buffered saline (DPBS) (Corning, #21-031-CM) was added to the inter-well space to avoid high evaporation. Three of the four 24-well plates were cultured in S.NEST at 10, 25 and 50 sec/cycle mixing rate for 4 days. The one remaining plate was cultured in a humidified standard incubator as control. Viable cell density and viability were measured with an automated cell counter (CURIOSIS, #FACSCOPE) on day 3

and day 4.

Statistical analysis

The data points were expressed as the mean±standard deviation (SD). Statistical analyses were performed using two-way ANOVA followed by Tukey's multiple comparisons test. A p value of less than 0.05 was considered a statistically significant difference. The P values are indicated listed as the following: P ≤ 0.05 (*), P ≤ 0.01 (**), P ≤ 0.001 (***) and P ≤ 0.0001 (****).

Results and discussion

In our study, CHO-S cells were seeded at various densities ranging from 2x10⁴ to 1x10⁶ cells/ml within four 24-well plates, as detailed in the "Materials and Methods" section. Three 24-well plates were employed in mixing culture at 10, 25, and 50 sec/cycle mixing rates in the NEST incubator. One plate was cultured in a humidified standard incubator as a control.


On day 3, the cell growth profile of the CHO-S cells was quantified by measuring viable cell density (VCD), relative expansion fold, and viability (**Fig 1**). As illustrated in **Fig 1A**, mixing groups with seeding densities at or exceeding 1x10⁵ cells/ml had significantly elevated VCD relative to the static control. Notably, for seeding densities at or higher than 5x10⁵ cells/ml, the 10 sec/cycle mixing speed significantly accelerated CHO-S cell proliferation compared to slower mixing speeds. Moreover, analyzing the relative expansion fold (**Fig 1B**) revealed that for seeding densities at or greater than 5x10⁴ cells/ml, dynamic mixing yielded superior expansion folds than static conditions.

These observations are consistent with our prior findings, underscoring the advantages of the NEST patented mixing culture, which provides homogenous nutrient and cell distribution uniformity and increases dissolved oxygen levels for cells. Nonetheless, viability measurements (**Fig 1C**) indicated that when cells were seed at the highest initial seeding density tested (1x10⁶ cells/ml), a faster mixing speed of 10 sec/cycle might have the drawback of lowering viability.

Measurements on day 4 echoed the results and trends observed on day 3 in VCD, expansion fold, and viability (**Fig 2**). Nonetheless, by comparing

results from day 3 and day 4, we can observe that high mixing speed may be beneficial at lower seeding densities when the culture is extended in duration. Specifically, on day 3, the 10 sec/cycle mixing speed conditions outperformed the 50 sec/cycle conditions at seeding densities at or higher than 5×10^5 cells/ml. However, on day 4, this outperformance occurred for seeding densities at or higher than 2×10^5 cells/ml. However, caution is warranted: although rapid mixing amplified the expansion fold for seeding densities at or higher than 2×10^5 cells/ml on day 4, decreased viability was also observed at higher seeding densities and mixing rates.

An overview of the cell growth improvement across different seeding density and mixing speed conditions is shown in **figure 3** and compared against the static culture controls. The chart shows that a NEST dynamic mixing culture could significantly accelerate cell growth at seeding densities exceeding 5×10^4 cells/ml in a 24-well plate. Collectively, these findings demonstrates that beyond the NEST mixing speed, initial cell seeding density also profoundly impacts cell growth, viability, and culture duration, warranting careful calibration.








Figure 1: Cell growth in different initial seeding density and mixing speed conditions after 3 days of 24-well plate culture. (A) Viable cell density. (B) Relative cell expression fold. (C) Cell viability.



Day4 Viability_24-well (NEST)




Figure 2: Cell growth in different initial seeding density and mixing speed conditions after 4 days of 24-well plate culture. (A) Viable cell density. (B) Relative cell expression fold. (C) Cell viability.

Seeding density (cells/mL)		2x10 ⁴	5x10 ⁴	1x10 ⁵	2x10 ⁵	5x10 ⁵	1x10 ⁶
Day 3	50s	-	+	+	+	+	+
	25s	-	+	+	+	+	+
	10s	-	+	+	+	++	++ Viability ↓

Seeding density (cells/mL)		1x10 ⁵	2x10 ⁵	5x10 ⁵	1x10 ⁶
Day 4	50s	+	+	+	+
	25s	+	+	+	+
	10s	+	++	++	++ Viability ↓

Figure 3: Overview of cell growth in different initial seeding density and mixing rate conditions in 24-well mixing culture. [-] no benefit over the static culture. [+] significantly improved cell expansion relative to static culture. [++] significantly improved cell expansion relative to mixing culture at lower mixing speed.

Conclusion

CHO cell culture optimization is pivotal in biopharmaceutical workflow. Utilizing the NEST microreactor, we explored the influence of initial cell seeding density and dynamic mixing speed on cell growth kinetics and viability in 24-well cell culture. Our data shows that high seeding densities, combined with high mixing speeds, enabled by the patented NEST microreactors, could provide optimal viable cell density (VCD). This observation is relevant to scaling bioprocesses where maximizing yield is imperative.

References

1. Wurm, F. M. (2004). Production of recombinant protein therapeutics in cultivated mammalian cells. *Nature Biotechnology*, 22(11), 1393-1398.
2. Kim, J. Y., Kim, Y. G., & Lee, G. M. (2012). CHO cells in biotechnology for production of recombinant proteins: current state and further potential. *Applied Microbiology and Biotechnology*, 93(3), 917-930.
3. Ozturk, S. S., & Palsson, B. O. (1991). Effect of medium osmolarity on hybridoma growth, metabolism, and antibody production. *Biotechnology and Bioengineering*, 37(10), 989-993.
4. Pillay, P., Maharajh, D., Panchal, M., & Chetty, S. (2017). The effects of cell culture parameters on cell-seeding density and monoclonal antibody production in CHO cells. *Biotechnology Progress*, 33(5), 1327-1339.



©2023 BICO AB. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of BICO is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO may refer to the products or services offered by other companies by their brand name or company name solely for clarity and does not claim any rights to those third-party marks or names. BICO products may be covered by one or more patents. The use of products described herein is subject to BICO's terms and conditions of sale and such other terms that have been agreed to in writing between BICO and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of BICO products in practicing the methods set forth herein has not been validated by BICO, and such nonvalidated use is NOT COVERED BY BICO'S STANDARD WARRANTY, AND BICO HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner BICO's terms and conditions of sale for the instruments, consumables or software mentioned, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by BICO that it currently or will at any time in the future offer or in any way support any application set forth herein.

Edited version: October 2023 | CBS_PUB_NEST_Tech-Note-04_Digital