


S/C.NEST® | Cross-Contamination-Free Operation Validated in S/C.NEST Platform


Shih-Pei (Betty) Lin, PhD, Missie Huang, BS
CYTENA Bioprocess Solutions, Taipei City, Taiwan



Introduction

In microscale bioprocessing, multi-well formats such as 24- and 96-well plates are essential for enabling high-throughput experimentation. The S/C.NEST microbioreactor developed by CYTENA BPS is designed to support these formats while minimizing the risk of cross-contamination. Using a patented air-driven mixing mechanism, S/C.NEST creates vertical fluid motion within each well through controlled airflow, delivered via a dedicated chamber and lid system. This airflow, sourced from both CO₂ cylinders and the surrounding atmosphere, is passed

through a HEPA filter to maintain sterility—similar to standard incubators.



This mechanism helps to preserve the physical separation of each well, supporting a better

containment even under continuous or vigorous mixing. This method enables consistent, effective mixing while maintaining isolation between wells — a critical advantage for sensitive assays or experiments requiring high reagent purity. Still, validating the system's containment performance remains important, particularly in demanding applications.

This study was designed to assess well-to-well cross-contamination risk during active mixing in the S/C.NEST platform by utilizing a visible dye in a checkerboard pattern across a standard 24-well plate.

To detect even minimal cross-contamination, the experiment uses a series of diluted Trypan Blue standards. These standards provide a reference to quantify dye presence in wells and surrounding areas through optical density measurements, allowing sensitive detection of any unintended dye transfer during mixing.


Materials and methods

A. Plate Setup and Reagents

A standard 24-well tissue culture plate was prepared using a checkerboard pattern in which alternating wells were filled with either:

- 1,400 μL of double-distilled water (ddH_2O), or
- 1,400 μL of a 0.2% Trypan Blue solution (prepared by diluting 0.4% stock 1:1 in ddH_2O).

This layout ensured that every Trypan Blue-filled well was completely surrounded by water-only "test wells". In addition to assessing well-to-well contamination, a separate focus of this experiment was to evaluate the potential for dye transfer between the moat (space between the wells) and the wells.



B. Device Configuration

The plate was incubated in a C.NEST equipped with

an X.NEST Lid, and operated under the following conditions:

- Temperature: 37 $^{\circ}\text{C}$
- CO_2 concentration: 5%
- Humidity: 90%
- Mixing Mode: Continuous with suspension boost
- Mixing Interval: Every 10 seconds

C. Experimental Conditions


To assess the potential for cross-contamination, experiments were performed under two different operational conditions:

1. Frequent lid openings — to simulate real-world sampling with:
 - Manual lid handling by an experienced user
 - EZHold-assisted lifting to reduce lid contact and droplet formation.
2. Long-term culture — to evaluate whether the platform's in-well mixing action alone causes cross-contamination over extended periods.

EZHold is a mechanical lid support tool that features a knob-operated lift mechanism that minimizes direct contact between the lid and the microwell plate surface. This reduces the chance of droplet formation or lid-edge residue that may result in cross-contamination, especially during repeated sampling or component additions.

In this experiment, lid handling was performed by an experienced S/C.NEST user. While the results presented here reflect a controlled scenario with low risk of contamination, it is important to note that manual operation by less experienced users may introduce higher variability.

To simulate real-world sampling conditions, the X.NEST lid was repeatedly opened and closed at defined intervals over a 3-hour culture period. At each sampling point, 200 μL was withdrawn from every well and transferred to a 96-well plate for measurement of optical density at 595 nm (OD_{595}) using a spectrophotometer. Trypan Blue served as the visual tracer for potential cross-contamination, and a standard curve was prepared in advance to quantify even low-level dye transfer.



D. Sensitivity and Detection Threshold

To accurately quantify potential cross-contamination, a series of Trypan Blue standards with known dilutions was prepared prior to the experiment. Optical density at 595 nm (OD_{595}) was measured for each standard. This curve served as a reference to estimate dye concentrations in samples taken from test wells and the surrounding moat during the experiment.

To ensure accurate OD measurements, each well was sampled twice—200 μ L was withdrawn twice as two separate replicates. Additionally, each 96-well plate was scanned three times to confirm consistency across readings and reduce variability.

By comparing OD_{595} readings of test samples to the standard curve, we could detect and quantify minimal amounts of dye — as little as a fraction of a microliter.

Results and discussion

Serial Dilution Curve

A standard dilution series of Trypan Blue was established before each experiment to correlate OD_{595} values with known dye concentrations. At each dilution, OD_{595} was measured for both the dye and a ddH₂O control blank. The blank's OD_{595} value was subtracted from the known dilution's OD_{595} values to obtain the corrected dye absorbance.




Figure 1: Trypan blue dilution curve

Although multiple dilution ratios were evaluated, this paper highlights the 100,000 \times , 10,000 \times , and 1,000 \times dilutions for clarity and relevance.

The expected amount of 0.4% Trypan Blue in each case is:

- At 100,000 \times dilution, up to 0.014 μ L of 0.4% Trypan Blue may be present in the solution.
- At 10,000 \times dilution, up to 0.14 μ L is present.
- At 1,000 \times dilution, up to 1.4 μ L is present.

To put this into perspective, a single water droplet from a pipette tip or condensation typically measures around 1–2 μ L. This volume is equal to or greater than the dye content in the 1,000 \times dilution, and substantially higher than in 10,000 \times or 100,000 \times . Therefore, even a minor cross-contamination event — such as the transfer of a single droplet of 0.2% Trypan Blue — would be detectable using this dilution curve. This makes it a sensitive and reliable reference for contamination monitoring.

Frequent Lid Openings – Manual Handling

The experiment began with 2 hours of continuous mixing, followed by a sequence of lid openings with 15-minute mixing intervals. In total, four lid openings were performed: the first after 2 hours, and the remaining three at 15-minute intervals after.

After each lid opening, samples were transferred from the test wells into a 96-well plate for OD_{595} measurement. The OD value of the ddH₂O blank (control) was subtracted from each measurement to yield a corrected OD value. The results obtained after the fourth lid opening is shown below:

OD value					
	0.037		0.040		0.037
0.039		0.038		0.038	
	0.037		0.036		0.035
0.036		0.036		0.035	
Moat:					0.037
(OD value) - (blank OD value)					
	0.001		0.004		0.001
0.003		0.002		0.002	
	0.001		0		-0.001
0		0		-0.001	
Moat:					0.001

Table 1: Test wells corrected OD_{595} value - Frequent opening using manual handling

Across all time points and lid openings, corrected OD₅₉₅ values in the ddH₂O test wells remained low—ranging between 0.000 and 0.006—with no progressive increase observed over time.

These results were compared with the Trypan Blue standard dilution curve:

- At 100,000x dilution, the OD₅₉₅ was 0.000.
- At 10,000x dilution, the OD₅₉₅ was 0.007.
- At 1,000x dilution, the OD₅₉₅ was 0.0825.

Since the test well readings remained below 0.006, any potential dye presence was below the detection limit defined by the 10,000x dilution, confirming the absence of meaningful cross-contamination.

Frequent Lid Openings – EZHold Handling

To assess whether EZHold impacts the risk of cross-contamination during frequent lid openings, a parallel test was conducted using the same dye-based approach. The experiment began with 1 hour of continuous mixing, followed by four sequential lid openings, each spaced by a 15-minute mixing interval.

OD value					
	0.037		0.036		0.036
0.036		0.036		0.036	
	0.036		0.036		0.036
0.036		0.036		0.036	
Moat:					0.036

(OD value) - (blank OD value)					
	0.001		0.000		0.000
0.000		0.000		0.000	
	0.000		0.000		0.000
0.000		0.000		0.000	
Moat:					0.000

Table 2 : Test wells corrected OD₅₉₅ value - Frequent opening using EZHold

OD₅₉₅ measurements of the ddH₂O test wells remained consistently low throughout the experiment, ranging between 0.000 and 0.001, indicating no measurable dye transfer.

Similar with the manual handling test, comparisons were made with the Trypan Blue dilution curve to evaluate detection sensitivity:

- At 100,000x dilution, the OD₅₉₅ was -0.0005.
- At 10,000x dilution, the OD₅₉₅ was 0.007.
- At 1,000x dilution, the OD₅₉₅ was 0.0825.

Since the test well readings never exceeded 0.001, the dye concentration remained well below the detection threshold defined by the 1:20,000 dilution.

Long-Term Culture – EZHold Handling

To simulate extended culture conditions, three separate plates were incubated with continuous in-well mixing for different durations: one for 3 days, one for 7 days, and one for 14 days. At the end of each respective culture period, the lid was opened, and the wells were assessed for potential dye transfer. In addition to the dye dilution curve, a static (no-mixing) control group was included to serve as a baseline for comparison.

OD value					
	0.037		0.036		0.036
0.036		0.037		0.036	
	0.037		0.036		0.036
0.036		0.036		0.036	
Moat:					0.036

(OD value) - (blank OD value)					
	0.00125		0.000		0.000
0.001		0.002		0.001	
	0.001		0.001		0.001
0.001		0.000		0.001	
Moat:					0.001

Table 3 : Test wells corrected OD₅₉₅ value - Long-term static culture

OD value					
	0.036		0.036		0.037
0.037		0.036		0.036	
	0.035		0.035		0.035
0.037		0.035		0.036	
Moat:					0.036

(value) -					
	0.0001		0.001		0.002
0.001		0.001		0.001	
	0.000		0.000		0.000
0.001		0.000		0.001	
Moat:					0.001

Table 4 : Test wells corrected OD₅₉₅ value - Long-term mixing culture

OD₅₉₅ measurements from the ddH₂O test wells across all plates—regardless of duration—remained low, ranging from 0.000 to 0.002, suggesting no evidence of dye cross-contamination even after two full weeks of continuous mixing.

Detection sensitivity was again benchmarked against the Trypan Blue dilution curve:

- At 100,000x dilution, the OD₅₉₅ was -0.0005.
- At 10,000x dilution, the OD₅₉₅ was 0.0085.

- At 1,000x dilution, the OD₅₉₅ was 0.096.

Even under long-term mixing conditions and with periodic lid opening, OD values in the test wells remained far below the lowest detectable levels on the curve.

Conclusion

This study evaluated whether cross-contamination occurs in the C/S.NEST platform under typical and laboratory conditions, using a Trypan Blue dilution curve as a sensitive detection reference. The dilution curve confirmed that even extremely small volumes of dye—equivalent to a single droplet—can be reliably detected using OD₅₉₅ measurements, providing a sensitive benchmark for monitoring contamination.

Throughout all tests, including manual and EZHold-assisted lid openings as well as long-term culture up to 14 days with continuous mixing, OD₅₉₅ values from ddH₂O test wells remained consistently below the detection threshold. No accumulation over time was observed, and values stayed lower than those associated with the 10,000× dilution (0.14 µL dye in 1400 µL), indicating no meaningful dye transfer.

Importantly, all OD₅₉₅ measurements taken from the moat—the space between wells—also remained at or near zero, further confirming that no cross-contamination occurred from the moat into neighboring culture wells when the platform was handled carefully.

Together, these results validate that the C/S.NEST platform maintains well isolation and spatial separation, even during frequent lid access and extended use. This makes it a robust option for contamination-sensitive applications in cell culture and beyond.



©2023 BICO AB. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of BICO is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO may refer to the products or services offered by other companies by their brand name or company name solely for clarity and does not claim any rights to those third-party marks or names. BICO products may be covered by one or more patents. The use of products described herein is subject to BICO's terms and conditions of sale and such other terms that have been agreed to in writing between BICO and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of BICO products in practicing the methods set forth herein has not been validated by BICO, and such nonvalidated use is NOT COVERED BY BICO'S STANDARD WARRANTY, AND BICO HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner BICO's terms and conditions of sale for the instruments, consumables or software mentioned, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by BICO that it currently or will at any time in the future offer or in any way support any application set forth herein.

Edited version: July 2025 | CBS_PUB_SNEST_Tech-Note-03_Digital